


CHROMABOND® PFAS



New solution for the enrichment of per- and polyfluoroalkyl substances

Dr. Max Bielitz, 11. November 2020

Agenda




Introduction



**Current solutions and
limitations**



New product solution



**Potential customers/markets/
competitors**



Summary



Introduction

Properties – daily use – guidelines



Introduction



Properties

- ~ 4730 compounds known according to OECD*
- All are nonnatural but man-made (used since the 1940s)
- Resistant against aggressive chemicals
- Mostly non-biodegradable (persistent) and accumulate in the environment and the food chain
- Often toxic, some can cause cancer, some are officially banned



* OECD: Toward a New Comprehensive Global Database of Per- and Polyfluoroalkyl Substances (PFASs): Summary Report on Updating the OECD 2007 List of Per- and Polyfluoroalkyl Substances (PFASs), Series on Risk Management, No. 39, ENV/JM/MONO(2018)7.



Hazardous compounds


Introduction




Properties and examples of PFAS

- Organic compounds
 - Carbon chain: hydrogen is substituted by fluorine
 - Carbon-fluorine bond very strong
 - Structure: hydrophobic, lipophobic chain + hydrophilic „head“
- tensid-like → water-, dirt- and fat-repellent (non-sticky)

Perfluorosulfonic acid



Perfluorooctanoic acid




Introduction



Appearance/daily use


- Fire-fighting foam
- Fiber coating
- Textile coating, e.g. seat covers, carpets, outdoor clothing
- Cookware
- Paper finishing
- Food packaging, e.g. pizza cartons, paper cups
- Building material, e.g. water resistant lacquer



Abundant in the (global) environment


Distribution of PFAS in the environment-examples

Precursor




Fire fighting with foam
Per/polyfluorinated compound
in fire fighting foam (etc.)

Dead-end products



$n = 3$ PFBA, Perfluorobutyric acid
 $n = 4$ PFPeA, Perfluoro pentanoic acid
 $n = 5$ PFHxA, Perfluoro hexanoic acid





Current solutions and limitations





Current solutions and limitations

Guidelines

- EPA 533, EPA 537 (**2008**), EPA 537.1 (**2020**)
 - PFAS from water
- DIN 38407-42 or DIN 38414-14 (**2011**)
 - 10 PFAS from water, soil



For now few official methods were developed



Current solutions and limitations

Weak anion exchanger – CHROMABOND® HR-XAW

- Recommended by DIN 38407-42 and EPA 533
- Methods are often sufficient for short-chain PFAS with acid functions
- But not for other PFAS substance classes
 - Alcohols
 - Alkyl ethers
 - Sulfonamides etc.

CHROMABOND® HR-XAW

Technical data

Weak anion exchanger based on polystyrene-divinylbenzene copolymer (PS/DVB)

SPE mode: Ion exchange and reversed phase (mixed-mode)

Interactions: Ionic, hydrophobic and π - π

Particle shape: Spherical

pH stability: 1–14


Particle size: 85 μ m and 45 μ m

Pore size: 55–65 Å

Specific surface: 850 m²/g

RP capacity: 350 mg/g (caffeine in water)

Exchange capacity: > 0.5 meq/g, pKa = 6



Recommended application

- Perfluorinated surfactants
- Acidic compounds like sulfonates
- Active ingredients from heavily matrix-contaminated samples, e. g., urine, plasma, serum
- Strong acids with pKa < 1

Good to know


A possible replacement for:

- Oasis® WAX
- Strata™ XAW



Not sufficient for all PFAS substance classes

Per- and polyfluoroalkyl substances (PFASs)



Per- and polyfluoroalkyl substances (PFAs) Source:

Buck, R.C.; Franklin, J.; Berger, U.; Conder, J.M.; Cousins, I.T; Voogt, P.d.; Jensen, A.A.; Kannan, K.; Mabury, S.A; Leeuwen, S.P.v (2011) Perfluoroalkyl and Polyfluoroalkyl Substances in the Environment: Terminology, Classification, and Origins. Integr.Environ.Assess.Manag Vol. 7 No. 4 pg. 513-541.

OECD, 2013, OECD/UNEP Global PFC Group, Synthesis paper on per- and polyfluorinated chemicals (PFCs),

Environment, Health and Safety, Environment Directorate, OECD. (http://www.oecd.org/env/ehs/riskmanagement/PFC_FINAL-Web.pdf, retrieved: 01 Jul 2020




Current solutions and limitations

Increasing interest

Continuous research on:

- Health effects
- Occurance
- Metabolism
- Reactant and Reagent
- **Analytical methods**

References (PFAS) in Scifinder



Scifinderⁿ : data retrieved: 13th Jul 2020



Need for advanced analytical methods and product solutions for monitoring



New product solution

Product and application note



New product solution



Special phase - CHROMABOND® PFAS

- Polymer-based combination phase with weak anion exchange functionality
 - Combines the advantages of different sorbents by using several retention mechanisms (dipole-dipole, ionic, hydrophobic, H-bond)
- Advantages
 - Solution for various PFAS substance classes
 - 28 PFAS can be enriched
 - Sorbent retention mechanisms according to DIN 38407-42, EPA 537.1 and 533 guidelines
 - High capacity
 - High recovery rates



New product solution



- One type of column:
 - 6 mL PP column
 - 300 mg adsorbent weight
- Standard packaging size: pack of 30
→ **REF 730283**
- BIGpack: pack of 250 (5 x 50 columns)
→ **REF 730283.250**
- Samples available
→ **REF 730283.Muster**



New product solution



CHROMABOND® PFAS – application note

- PFAS from water



New product solution



Chromatographic conditions

- Solid phase extraction
 - Column: CHROMABOND® PFAS, 6 mL, 300 mg
 - Conditioning: 10 mL 0.1% NH₃ in methanol, 10 mL methanol, 10 mL water
 - Sample application: 150 mL water sample with a flow rate of 5 mL/min
 - Washing (optional): 5 mL of 25 mM (NH₄)⁺ acetate buffer (pH 4.0) water with flow rate of 3 mL/min
 - Drying: 1 min with vacuum
 - Elution: 10 mL 0.1% NH₃ in methanol
 - Eluent exchange: Evaporate eluate to dryness at 40 °C under a stream of nitrogen and dissolve residue in 0.5 mL water / methanol (40:60, v/v)
- Subsequent analysis: LC-MS/MS with NUCLEOSHELL® RP 18plus (2.7 µm)

New product solution



Recovery rates

- For most PFAS: **80 – 100 %**
- Lower for
 - FOSA (sulfonamide) lower due to high polarity/low ionic interaction
 - PFTTrDA, PFTeDA, PFDoA due to solubility issues and adsorption effects on surfaces



Good and reproducible recovery rates for most analytes



Potential customers/markets/competitors



Potential customers/markets/competitors



- Contract laboratories (environmental)
- Textile and paper industry
- Chemical industry, producers of:
 - Fire extinguisher foams
 - Lubricants, emulsifying agents
 - Impregnation agents, etc.
- Galvanic industry
- Aviation
- Photo and semi-conductor industry
- ...






Summary




Summary



- PFASs are persistent, tensid-like, harzardous compounds
- Commonly used because of their non-sticky properties
- Introduced by textiles, fire extinguisher foams, paper cups etc.
- Used since the 1940s, but long-time neglected
- Found in the environment (water, foods, soil, animals and humans)
- **Current analytical methods for monitoring are limited**



Perfluorosulfonic acid



Perfluorooctanoic acid



Summary



- **CHROMABOND® PFAS:**

- Solution for the enrichment of a broad range of different PFAS (32) from water
- Several sorbent retention mechanisms according to DIN 38407-42, EPA 537.1 and 533 guidelines
- High reproducibility and recovery rates
- Application note available



- Samples available
- Customers: contract laboratories (environmental); textile, paper and chemical industry, etc.

Thank you for your attention!

Dr. Max Bielitzka | mbielitzka@mn-net.com | +49-2421-969-375

Image credits

- © Aaron Amat (2-6), vipman4 (2, 13-18), djama (2, 19-21), ad_stock (2, 22-24), nasir1164 (4), New Africa (6), MandicJovan (6), bymandesigns (6), Karyna Chekaryova (6), kegfire (6), vegefox.com (9), Coloures-Pic (14), BRAD (16), Olivier Le Moal (18), Gresei (20), nonameman (23)